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LEONE IMPLANT SYSTEM - INSTRUCTIONS FOR USE

Warnings and restrictions. The herewith furnished information is NOT intended for the patient, but exclusively for the medical professional, who will have to evaluate the contents in reference to every specific case to be treated, in order to furnish suitable information to the patient about the benefits and eventual side effects or complications. Leone Implant System must only be used by Dentists or Dental Surgeons experienced in implant techniques as well as in the related diagnostic methods, pre-operative planning and manufacture of prosthesis on implants. The hereunder instructions for use are finalized to underline the risks, the contraindications and the collateral effects of general character (a reference to all the possible circumstances of use is not possible) and well known at the moment of the placing of the device on the market. Nevertheless, the Leone Implant System requires the knowledge of the surgical and prosthetic procedures illustrated in the Leone Guidelines: Surgical Procedure and Prosthetic Procedure available either on Leone Implantology Catalogue or on the web site www.leone.it, which specify in detail the indications and the product restrictions of use. Considering that errors and deficiencies in patient evaluation, preoperative diagnosis, in the planning and the carrying out of the treatment may lead to implant failure, before use, it is strongly recommended to attend specific training courses with the purpose to reach a high degree of propaedeutical knowledge and practice in the use. It is not ascribed to Leone the duty of monitoring the procedures of use of the products for the implant-prosthetic treatment; such procedures must integrally be reverted under the responsibility of the clinician. Leone takes no responsibility for the application of the device and its workmanship nor for possible surgical or prosthetic misuse of the device even neither for failure, adverse reactions or damages occurred to the patient or to the dentist following the application of the product. In any case the use of exclusively original components, accessories and instruments in conformity to the instructions is strongly recommended. **U.S.A. warning.** "Caution: Federal law restricts this device to sale by or on the order of a physician." **Description of the product.** The dental implants, the healing caps and the low self-locking caps of the Leone Implant System are medical devices made of titanium for surgical use (ISO 5832-3). The definitive abutments are medical devices made of titanium for surgical use or gold alloy; the cover caps and the temporary abutments are made of polymers for implant use. The implants and the abutments are provided with specific accessories which include components and instruments for surgery, prosthetic dentistry, and dental laboratory. The accessories and the instruments can be made of titanium, stainless steel, and biocompatible plastic material. The implants, the abutments, and related accessories are available in different sizes in order to match individual anatomical variability. **Intended use.** The Leone Implant System is indicated for single or multiple tooth replacement, for use in terminal or intermediate edentulous sites in the mandible and/or the maxilla, and for totally edentulous arches. The system is designed to be surgically inserted in the bone structure of the mouth in order to replace missing teeth. It can work as a support for partial/total prosthetic restorations, either provisional or definitive, or as an anchorage system for removable prosthesis. **Contraindications.** All the pathological conditions, both local and systemic, which contraindicate an intervention of oral surgery, have to be considered contraindications valid for the application of dental implants as well. As information only, and provided that a specific clinical evaluation of the patient is strictly recommended, it is understood that the final choice to proceed or not with the treatment is an exclusive competence of the sanitary clinician, Leone implants should not be used in patients with recent myocardial infarction, haematological disorders, metabolic bone disorders, neoplastic pathologies diagnosed during the last 5 years, immunodepressed patients, previous bone radiotherapy or chemotherapy, patients with cardiovascular or metabolic problems which are not compensated, dismorphic patients, severe mental illness, xerostomia, diabetes or hyperglycaemia, transplanted patients, renal insufficiency, fibrous dysplasia, pharmacological treatments with anticoagulants, anticonvulsants, immunosuppressants, steroids, acute or chronic infectious illnesses, pregnancy, nursing, insufficient healthy bone volume and impossibility or unavailability of the patient for the augmentation technique, skeletal growth still in action, scarce oral hygiene, parodontitis, candidosis and pathologies of the oral mucosa, allergic reaction to titanium, insufficient arch space, occlusion and/or articulation defects, para-functional habits (clenching, bruxism, etc.), poor patient's compliance and motivation, abuse of drugs, alcohol and tobacco, acute, sub-acute or chronic maxillary osteitis, radicular residuals. It is understood that the final choice to proceed or not with the treatment is an exclusive competence of the sanitary clinician. The implants inserted in the maxillary arch should not pierce the maxillary sinus. The implants inserted in the mandibular arch should not touch, nor compress, or cut off the mandibular nerve. **Side effects.** Apart from the features of the utilized implant, the surgical application can be followed by edema, pain, haematomas, inflammations, limitation or alteration of the sensibility and oral functions, allergic reactions. Generally such effects are temporary and the patient has to receive a complete information about these side effects and about the necessity to give immediate communication to the physician when they would occur. **Complications and interactions.** Referable complications related to the Leone implant system are not known neither predictable. The surgical procedure, particularly if defectively performed, may cause hemorrhagic diseases, permanent paresthesia and/or dysesthesia, bone re-absorption at mandibular/maxillary crest level, infectious processes, adversely affected adjacent teeth, implant or prosthesis breakage, bone fractures, mucous ulcers, gingival hyperplasia, oroantral or oronasal fistulae, aesthetical problems, nervous lesions, inhalation or ingestion of parts or tools of the implant system. In the event that parts or fragments remain in the tissue after the breakage, or inhalation/ ingestion, the removal and the relative diagnostic checks are advised (keeping in mind that the devices are all radiopaque, with the only exception for the temporary abutments). **Packaging and sterility.** The Leone dental implants, the cover caps, the healing caps and the low self-locking caps are supplied gamma ray sterilized in glass containers. The product remains sterile if packaging is intact. When stored correctly, the sterility of the product is guaranteed up to the expiration date. Do not use the product after the expiration date printed on the label. Store the product in the original packaging, in a dry place, at room temperature and away from sun light. Leone takes no responsibility in case of re-sterilized dental implants and caps, regardless of the method of sterilization. Never use non-sterile or pre-used implants and caps. All abutments and prosthetic components are supplied non sterile and must be sterilized using a validated process, which should be constantly monitored and conforming to ISO 17665-1 "Sterilization of health care products - moist heat - Part 1: requirements for the development, validation and routine control of a sterilization process for medical devices". It is recommended to sterilize the components with steam autoclave using the following parameters: temperature: 121°C (250°F), pressure: 1 atm, minimum exposure time: 20 minutes. The sterilization packaging must be appropriate for the device and conforming to EN 868-5 "Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods". **Identification and traceability.** The dental implants, the cover and healing caps, the low self-locking caps, the abutments and the other prosthetic components are identified by means of the code and the lot number. Such information is essential to assure the traceability of the product. The package is provided with an adhesive label showing these data that is to be applied on the patient's record sheet. **Final precautions.** The patient should be informed about the precautions for the period after installation of the implant, during and at the completion of treatment, in order to prevent complications, and variations in the efficiency of the device. For instance, heavy physical activity, immediate post-surgical mechanical loading of the implant, and occlusal trauma should be avoided. A good level of oral hygiene should be achieved and maintained. Periodical check-ups should be performed.

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Below are the symbols used, where applicable, on the label of the products of the Leone Implant System.

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