


LEONEAMERICA

# ORTHO NEWS

A PERIODICAL FOR ORTHODONTIC PROFESSIONALS

SPRING/SUMMER 2019  
VOL 1 #41



**"LEAF**  
IS A GAME  
CHANGER"

**Yan Razdolsky**

DDS, BSD, LTD  
Diplomate of American  
Board of Orthodontists

**LEONEAMERICA**<sup>TM</sup>



SEE PAGE 5



MANUFACTURING  
ORTHODONTIC PRODUCTS  
SINCE 1934



# TABLE OF CONTENTS

- 5 Interview to Dr. Yan Razdolsky
  
- 11 Use of the Leaf Expander® in the treatment of transversal discrepancy in adults: a clinical case
  
- 19 Ceramic Bracket – Interactive SL (Self-Ligating)
  
- 20 Pediatric Orthodontics
  
- 21 Leaf Series Expander
  
- 22 Leone GoTo1 (Class II Corrector)
  
- 23 Where To Find Us: LeoneAmerica Collaborations

## Pedodontic Bands



Designed to allow for early treatment of patients with mixed or deciduous teeth.

## Self-Ligating Brackets InteractiveSL

**Inter. \ctiveSL**

D.B. Metal  
Brackets  
Self-Ligating




The specific shape of the clip and slot allow users to modulate the most appropriate level of friction force between bracket and wire.



**ORTHODONTICS  
FOR KIDS**

*“A CHILD  
WITHOUT SMILE  
IS LIKE A GARDEN  
WITHOUT FLOWERS”*

**PRODUCTS DEDICATED TO PEDIATRIC TREATMENT**



# Interview to Dr. Yan Razdolsky



*Dr. Yan Razdolsky, DDS, BSD, LTD, Orthodontist for Children and Adults Dr. Razdolsky is a Diplomate of the American Board of Orthodontics. His professional affiliations include the American and European Associations of Orthodontists, World Federation of Orthodontists, the Midwestern and Illinois State Societies of Orthodontists, the American Dental Society, the Illinois and Chicago Dental Societies, Academy of General Dentistry, and the Alpha Omega International Dental Fraternity. Dr. Razdolsky is also an active member of reputable Seattle Study club.*

Yan Razdolsky, long-term friend of Leone and LeoneAmerica, in this brief interview, shares his experience with the Leaf Expander, highlighting the clinical advantages but also the advantages introduced in the management and organization of his dental practice.

## **1. How many palatal expansion have you performed during your career?**

I have been in practice 32 years. We place great emphasis on early prevention. I would estimate over 40,000 RPEs to date.

## **2. How did you get to know about the Leaf Expander?**

I have met Gabriele Scommegna in mid - nineties. I was heavily involved in Distraction Osteogenesis, bone growth via slow stretching.

I owned 5 patents for mandibular, maxillary, and vertical alveolar distraction. Leone was the only company in the world who answered my call and have designed anterior activation screws for us. We have visited Leone facility in Firenze and were met with outmost hospitality! I have lectured on Distraction Osteogenesis all over the world but unfortunately the technique never gained wide spread popularity and my patents have long expired. Gabriele and I kept in touch throughout the years and became very good friends. So when we met at AAO Mid-winter meeting in Scottsdale 2018 he showed me Leone's new invention – Leaf Expander!

## **3. How many cases have you treated with the Leaf Expander?**

We use Leaf expander on most of our expansion cases now including SARPEs and Herbst. I would estimate over 100 Leafs in the first year of use. Enclosed a case study!

## **4. What was your main concern about the use of the Leaf Expander?**

I knew right away that this is a major breakthrough as far as the clinical efficiency is concerned. I could not wait to test the efficacy of Ni-Ti Leafs. Gabriele has shared with me clinical studies from Italy and I was confident the appliance will deliver. JCO article a bit later was also a confidence builder.

## **5. On average, how many visits does a patient need during a Leaf Expander treatment and how many during an RPE treatment?**

Let's consider the RPE appointments. If it's a Bonded RPE we would scan the case and download it on the laboratory server. With banded RPE we learned the hard way that it's much more precise to fit your own bands, take a rubber base impression and mail it to the lab. So this aspect does not vary between RPE and the Leaf. Delivery is no different either. However, with RPE you have



# Interview to Dr. Yan Razdolsky

to give patient few weeks to get used to, get them back to teach parents how to expand, and then continue to manage the case if parents do not perform or are not dexterous enough to deliver the results. We literary have had patients who were unable to turn the screw properly and ended up having patient come back on a weekly basis to activate 3-4 turns a week. So on the average using RPE vs Leaf can add 3-4 extra appointments and what about staff and doctors time?

## 6. Which are the main advantages of the Leaf Expander?

Improvement in clinical efficiency. You deliver the Leaf, see patient every 6-8 weeks and perform 15 turns every appointment. The whole process is streamlined and many times my staff is already done by the time I come to the chair.

No learning curve for parents or patients. You project much greater value for what they are paying. Why would you pay someone and have to perform the work yourself? People are way too busy now days. Using conventional RPEs adds xtra tasks for them.

No calls from the patients re: instructions, questions, difficulties, etc..

## 7. How did the Leaf Expander affect the management of your daily activities in the dental office?

To add to my summary in Q6., we see on the average 85 patients/day. If we can cut out few unnecessary appointments and save our patients an extra trip and aggravation it just makes the whole experience so much greater!

## 8. Are there any contraindications in the use of the Leaf (patient age, malocclusions class etc.)?

So far I have not found any. In SARPE we see patients every two weeks (instead of 6-8 weeks) because we have a narrow time frame to complete the process before the boney fusion. Patients love the fact that they do not have to activate it themselves as everything is so tender after the surgery. I would love to have a smaller Leaf to be used in tooth borne or TAD supported distallizers.

## 9. Did you have issues for the manufacturing and installation of the Leaf Expander?

We have discovered that during the soldering in the laboratory the plastic straps loosen up possibly due to heat. It makes the Leaf expand and therefore harder to deliver. So have been ligating the screw prior shipping them off to the lab along with the case.

## 10. Do you think that the use of the Leaf Expander needs peculiar expertise?

Very little learning Curve!

## 11. What about care givers' and patients overall satisfaction?

People loved us but they love us even more now! My staff is elated!

## 12. What would you tell to an orthodontist who is still suspicious about the use of the Leaf expander?

Once in a while an invention appears which fundamentally changes the whole game. Leaf is a game changer!

Start using it and you will quickly realize the benefit both to you and your patients!





*Gabriele Scommegna Research & Development Director Leone Spa,  
Dr. Razdolsky*




*Hyrax – Class III Ectopic #3*




*9 yrs old  
Class III Skeletal w/Unilateral crossbite*



*Hyrax – Class III Ectopic #3*



*Hyraax – Class III Ectopic #3*



*Leaf delivered May 2018 w/De-impactor spring and Facemask Hooks*





*3 Leaf Activations x 15 turns performed 6-8 weeks apart  
Facemask used at night*



*Debonded December 2018  
Complete resolution of crossbite, crowding., and underbite!*



*#3 Fully erupted*

# MEMORIA® LEAF SPRING ACTIVATED EXPANDER



Nickel Titanium MEMORIA®  
leaf springs allow the release of  
calibrated and continuous forces  
to promote the expansion of the  
maxillary arch.

**ok**  
ORTHODONTICS  
FOR KIDS

# Use of the Leaf Expander® in the treatment of transversal discrepancy in adults: a clinical case

Authors: Grecolini M.E., Casali A., Celli D. and Mele G.

Original title: *Utilizzo del Leaf Expander® nel trattamento del deficit trasversale mascellare dell'adulto: un caso clinico, Leone Bolletino number 101, June 2018.*

## Introduction

Rapid palatal expansion represents the gold standard procedure used for the correction of cross-bites. The ossification of the mid-palatal suture is a key factor in the decisional process for the implementation of the correct treatment plan.

Although some studies (1) report that the ossification of the mid-palatal suture can be observed between 15-19 years of age of the patient, some others state that it is not noticeable at 27, 32, 54 and 71 years (1, 2, 3).

Therefore, considering that the biological age is not a valid decisive factor (2,3,4), the vertebrae maturation index and the use of the CBCT (Cone Beam Computed Tomography), both in pubertal and post-pubertal age, are considered the most reliable methods (5,6,7,8).

If it is not possible to perform the mid-palatal suture opening (or only partially), dento-alveolar expansion represents the only achievable result using maxillary expanders. The dento-alveolar expansion is gained in the 39-49% of the cases and represents from the 6% to the 13% of the total expansion (10,11).

## Treatment Plan Summary

In light of the adult age of the patient and following the specific patient request to avoid surgical procedure and/or tooth extraction, the severe skeletal malocclusion has been treated with a no-surgical approach.

The patient had a severe mono-lateral cross-bite with dental crowding both in the upper and in the lower arch, misplacing of the inter incisive line, Class I occlusion at left side while Class II, both molar and canine, malocclusion on the other quadrant (Fig. 1).

A Leaf Expander 900 g has been used to treat the maxillary transversal discrepancy.

Leone Leaf Expander is a palatal expander that, in adults, can induce a modification of the transversal dimension mainly

through a dento-alveolar remodeling. This result is achievable thanks to the super-elastic properties of the Nickel Titanium unique leaf shaped springs, which release controlled and continuous force.

Since the patient had a asymmetrical transversal discrepancy, the use of cross elastic traction for the right side has been necessary. A lingual arch has been used for the control of the lower anchorage.

In order to guarantee an optimal control of the torque, Straight Wire appliance was used: upper laterals brackets have been placed in advanced stage of therapy, rotated of 180° and then normalized at the end of the therapy.



Fig. 1 Beginning of the treatment: X-Rays

---

# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

---

## Clinical facial analysis (Fig. 2)

- Facial asymmetry with skeletal high angle, prominent mandibular body highlighted by the slimness of the patient,
- Irregular smile with light labial incompetence,
- Noticeable altered orbicularis musculature.

**Functional analysis** – The patient refers recurrent muscular-tensile pain. Presence of sporadic TMJ bilateral click and pain. The misalignment of the midlines is not noticeable when the mouth is open, highlighting therefore a mandibular shift, probably due to an initial incoordination of the condyle-meniscus. As per patient decision, no further MRI investigation was performed.

**Intra-oral analysis (Fig. 3)** – Mono lateral cross bite (right), upper and lower severe dental crowding. Misalignment of the mid-line, asymmetric molar and canine occlusion: Class I at left and Class II at right. Cross palatal position of the upper laterals. No OVB, No OVJ.



Fig. 2 Beginning of the treatment: extra-oral pictures



# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

## Model Analysis (Fig. 4)

### Upper arch:

- Absence of third molars
- Asymmetric arch
- Severe dental crowding
- Palatal position of upper laterals and vestibular ectopic canines
- Both side molars rotated.

### Lower Arch

- Absence of third molars
- Severe asymmetric arch
- Severe grade of dental crowding
- vestibular ectopic position of left canine
- Rotations



Fig. 3 Beginning of the treatment: Intra-oral pictures




Fig. 4 Cast models: before and after treatment

# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

## Occlusal relationships-sagittal plane

- Molar and canine occlusion: Class I at left and Class II at right
- No OVJ
- Deep Spee Curve

## Occlusal relationships-frontal plane

- No OVB,

## Occlusal relationships-transversal plane

- Severe reductions of maxillary and mandibular transversal diameter
- Cross bite at upper laterals due the palatal positions

**Treatment Plan** — Skeletal age of the patient and the malocclusion type were crucial to determine the treatment plan. Due to the severe reduction of the maxillary transversal dimension and the lateral cross-bite, an important dento-alveolar remodeling was necessary and moreover it was necessary to torque heavily the left bicuspid and first molar, and lightly the right ones.

A Leone Leaf Expander has been used that release 900 gr and allows 6mm of expansion (A2704-06) (Fig 4.a).

The screw activation has been managed fully by the doctor in three treatment sessions during the whole cycle of treatment (Tab.1).

Thanks to this protocol, it was possible to apply constant pre-determined and calibrated force that allowed to fully control the expansion movement and the vestibular inclination of "pillar" teeth.



Fig. 4.a: Leone Leaf Expander 900g, 6mm (ref. A2704-06)

Tab. 1: protocol for Leaf Expander activation

EXPANDER ACTIVATION GUIDELINES	LEAF expander cementation	FIRST visit	SECOND visit	THIRD visit
<b>6 mm</b> <b>A2703-06 - A2704-06</b>	spring activation 0 turns	after 6 weeks 10 turns	after 4 weeks 10 turns	after 4 weeks 10 turns
<b>9 mm</b> <b>A2703-09 - A2704-09</b>	spring activation 0 turns	after 8 weeks 15 turns	after 6 weeks 15 turns	after 6 weeks 15 turns

# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

As reported in the literature, assuming the device is accurately manufactured with an optimal fit on the “pillar” teeth a high control of vestibular tooth inclination, by means of a corporeal movement in the vestibular direction, could be achieved (13,14,15,16).

Considering the severe maxilla asymmetry, the use of cross bite elastics, on the right, side has been necessary. The excellent compliance of the patient allowed to obtain a good correction in a relatively short time.

Upper laterals were corrected using sequential bonding focusing on the correct positioning of the roots. In particular, the two brackets were placed upside-down in order to exploit a wider correction of the root vestibular torque.

The lower arch was challenging due to the severe dental crowding and the cuspid ectopic vestibular position. The correction has been performed initially keeping a high control of the anchorage and paying attention to the anterior teeth pro-inclination.

Vertical and asymmetric elastics were used to obtain the correction of the malocclusion class and the correct OVB. The pink esthetics and the symmetry of the marginal gums has been obtained with the re-contouring of the gingival margins using a laser diods (Fig. 5).

## Results

The whole treatment took 1 year and 10 months. The results can be summarized as follow.

**Skeletal** — For what concerns the cephalometric parameters, the numerical changes are noticeable in the improvement both of the malocclusion class with a light reduction of the AN/Pg due primary to the mandibular advancement and the skeletal high angle with the reduction of the SN/GoGn angle (Tab.2).

The correction of the maxillary transversal dimension has been obtained thanks to a dento-alveolar remodeling using a Leaf Expander. Since the therapy has been performed on a not growing adult patient, the skeletal changes are minimal (Fig. 6).

**Soft tissues** — The facial features are noticeably changed with an improvement of the facial symmetry and smile. The musculature is relieved and the labial competency enhanced.

**Dentals** — A good dental alignment has been obtained with Class I occlusion both at the canine and molar level. The mid lines are aligned, OVJ and OVB are within normal value.

After the removal of the orthodontic appliance, the patient underwent to a professional dental whitening procedure and in order to improve the smile, a laser gingivectomy of the upper arch has been performed.

The smile present a valid esthetic line and is more expressive and harmonious with the patient face (Fig.7, Fig.8, Fig.9 and 10).




Fig. 5 Esthetic re-contouring of the gingival borders and elastics



# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

Tab. 2: Summary table of dental and cephalometric parameters

OCCLUSAL RELATIONSHIPS-SAGITTAL PLANE			
Maxilla position S.N/A	82° +/- 3,5°	76°	76,5°
Mandibular position S.N/ Pg	80° +/- 3,5°	71°	73°
Sagittal Inter-Maxillary relations A.N/Pg	2° +/- 2,5°	5°	3,5°
OCCLUSAL RELATIONSHIPS-FRONTALE PLANE			
Palatal plane inclination S.N/ANS.PNS	8° +/- 3,0°	16,5°	14°
Mandibular plane angle S.N/Go.Gn	33° +/- 2,5°	44,5°	42,5°
Vertical Inter-Maxillary relation ANS.PNS/Go. Gn	25° +/- 6,0°	28°	28,5°
DENTAL BASAL RELATIONSHIP			
Upper Incisor inclination +1/ANS.PNS	110° +/- 6,0°	113°	111°
Lower Incisor inclination -1/Go.Gn	94° +/- 7,0°	102°	95°
Lower Incisor compensation -1/A.Pg. (mm)	2 +/- 2 mm	5 mm	XXX
DENTAL RELATIONS			
Overjet (mm)	3,5 +/- 2,5 mm	0	XXX
Overbite (mm)	2,5 +/- 2,5 mm	0	2 mm
Inter incisive angle	132° +/- 6,0°	117°	124,5°





Fig. 6: Comparison of X-Rays Images and cephalometric tracing before and after treatment.



---

# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case

---



*Fig. 7 End of the treatment: extra-oral pictures*



*Fig. 8 End of the treatment: Intra-oral pictures*

# Use of the Leaf Expander<sup>®</sup> in the treatment of transversal discrepancy in adults: a clinical case




Fig. 9 Comparison before-after treatment digital upper arch cast




Fig. 10 Digital cast palatal rugae superimposition

## Reference

1. Persson M, Thilander B, Palatal suture closure in man from 15 to 35 years of age. *Am J Ortho.* 1977;72:42-52.
2. Knoup B, Yildizlan F, Wehrbein H, Age -related changes in the mid palatal suture. *J. Ortofac. Orto.*2004; 467-74.
3. Korbmacher H, Schilling A., Puschel K., Amling M., Kahl-NieKe B., Age dependent three dimensional micro computer tomography analysis of the human mid palatal suture. *J. Orofac. Orthop.* 2007; 68:364-76.
4. Persson M, Magnusson B.C, Thilander, B. Sutural closure in rabbit and man: morphological and histochemical study(Article). *Journal of Anatomy* Volume 125, Issue 2, 1978.
5. Baccetti T, Franchi L, Mc Namara J.A., The cervical vertebral maturation (CVM) method for the assessment of optimal treatment timing in dentofacial orthopedics. *Seminar Orthod.* 2005; 11: 119-29.
6. Franchi L, Baccetti T, Mc Namara J.A. Jr, Mandibular growth as related to cervical vertebral maturation and body Height. *Am. J. Orthop.* 2000; 118:335-40.
7. Perinetti G, Caprioglio A, Contardo L., Visual assessment of cervical vertebral maturation stages: a study of diagnostic accuracy and repeatability. *Angle Orthod.* 2014; 84:951-8.
8. Angelieri F, Cevidanes L.H., Franchi L, Goncalves J.R., Benavides E., Mc Namara J.A. Jr, Midpalatal suture maturation: classification method for individual assessment before rapid maxillary expansion. *Am J. Orthod. Dentof. Orthop.* 2013; 144:759-69.
9. Angelieri F, Franchi L, Cevidanes L.H.S., McNamara J.A., Diagnostic performance of skeletal maturity for the assessment of midpalatal suture maturation. *Am J Orthod Dentofacial Orthop* 2015;148:1010-6.
10. Garrett B.J., Caruso J.M., Rungcharassaeng K., Farrage JR., Kim Js, Taylor GD, Skeletal effects to the maxilla after rapid maxillary expansion assessed with cone-beam computed tomography. *Am J Orthod Dentofacial Orthop* 2008; 134:8-9.
11. Weissheimer A, deMenezes L.M., Mezomo M, Dias D.M., deLima E.M., Rizzato SM. Immediate effects of rapid maxillary expansion with Haas-type and hyrax-type expanders: a randomized clinical trial. *Am J Orthod Dentofacial Orthop* 2011;140:366-76.
12. Grünheid T, Larson C.E., and Larson B.E., Midpalatal suture density ratio: A novel predictor of skeletal response to rapid maxillary expansion. *Am J Orthod Dentofacial Orthop* 2017;151:267-76.
13. Lanteri C., Lerda F., Francolini F. L'espansore lento ammortizzato (E.L.A.): Un nuovo apparecchio di espansione mascellare *Boll. Inform. Orto.* 4: 22-28 2005.
14. Lanteri C., Beretta M., Lanteri V. L'espansore lento ammortizzato (E.L.A.) *Boll. Inform. Orto.* 79:11-20 2007.
15. Lanteri C., Beretta M., Lanteri V. L'utilizzo dell'E.L.A. nell'espansione mascellare *Dent. Trib. III*, 7:6-12 2007.
16. Lanteri C., Lanteri V., Gianolio A., Beretta M., Cherchi C., Franchi L. A new way for no compliance palatal expansion: the Leaf Expander *JCO* 90: 552-560 2016.



# D.B. METAL BRACKETS




## Interactive control

The specific shape of clip and slot of **InterActive** brackets allows the user to modulate the most appropriate level of friction force between bracket and wire, depending on the needs of the various stages of treatment.



## Passive phase

First stage round arches are not bound by the clip: the low friction will facilitate the process of alignment and leveling.



## Interactive phase

Rectangular arch wires, used for space closure, rotation, and torque control, work to elastically deform the clip for the biomechanical control necessary in these stages of treatment.



## Active phase

Rectangular arches for finishing and detailing completely fill the slot while going into active contact with clip: this allows the exploitation of elastic properties in order to obtain minimal movements for finishing of treatment.

## Easy open/close

The clip has a central hole and does not require any special tool for opening and closing.



## Opening

Insert the tip of a probe or utility tool into the hole in the clip and move towards the occlusal plane.



## Closing

Slide the clip with a slight pressure towards the gums using a tool tip or even just a finger.

### InterActive Self Ligating D.B. Brackets ROTH System

<b>F1100-11</b>				+12°	+5°	1
<b>F1100-21</b>				+12°	+5°	1
<b>F1100-12</b>			2.8	+8°	+9°	2
<b>F1100-22</b>				+8°	+9°	2
<b>F1100-13</b>			2.6	-2°	+13°	3
<b>F1100-23</b>				-2°	+13°	3
<b>F1100-14</b>			2.8	-7°	0°	4
<b>F1100-24</b>				-7°	0°	4
<b>F1100-15</b>			2.8	-7°	0°	5
<b>F1100-25</b>				-7°	0°	5
<b>F1100-41</b>			2.8	0°	0°	1
<b>F1100-31</b>				0°	0°	1
<b>F1100-42</b>			2	0°	0°	2
<b>F1100-32</b>				0°	0°	2
<b>F1100-43</b>			2.5	-11°	+7°	3
<b>F1100-33</b>				-11°	+7°	3
<b>F1100-44</b>			2.8	-17°	0°	4
<b>F1100-34</b>				-17°	0°	4
<b>F1100-45</b>			2.8	-22°	0°	5
<b>F1100-35</b>				-22°	0°	5

Packs of 5

### Interactive Self Ligating D.B. Brackets Kit ROTH System

1 case 20 Brackets		10 kits of 1 case	
<b>F1100-91</b>	.022"x .030"	<b>F1101-91</b>	.022"x .030"



### InterActive Self Ligating D.B. Brackets MBT\* System

<b>F1102-11</b>				+17°	+4°	1
<b>F1102-21</b>				+17°	+4°	1
<b>F1102-12</b>			2.8	+10°	+8°	2
<b>F1102-22</b>				+10°	+8°	2
<b>F1102-13</b>			2.6	-7°	+8°	3
<b>F1102-23</b>				-7°	+8°	3
<b>F1102-14</b>			2.8	-7°	0°	4
<b>F1102-24</b>				-7°	0°	4
<b>F1102-15</b>			2.8	-7°	0°	5
<b>F1102-25</b>				-7°	0°	5
<b>F1102-41</b>			2.8	-6°	0°	1
<b>F1102-31</b>				-6°	0°	1
<b>F1102-42</b>			2	-6°	0°	2
<b>F1102-32</b>				-6°	0°	2
<b>F1102-43</b>			2.5	-6°	+3°	3
<b>F1102-33</b>				-6°	+3°	3
<b>F1102-44</b>			2.8	-12°	+2°	4
<b>F1102-34</b>				-12°	+2°	4
<b>F1102-45</b>			2.8	-17°	+2°	5
<b>F1102-35</b>				-17°	+2°	5

Packs of 5

### Interactive Self Ligating D.B. Brackets Kit MBT\* System

1 case 20 Brackets		10 kits of 1 case	
<b>F1100-92</b>	.022"x .030"	<b>F1101-92</b>	.022"x .030"




\*MBT is a 3M Unitek Trademark

All the orthodontic brackets illustrated in this brochure are not intended to be a duplication of any other existing bracket nor does Leone Spa imply that they are endorsed by the above mentioned doctors or Schools.



**ORTHODONTICS FOR KIDS (OK) PEDIATRICS**

Early treatment options are continuously a popular topic in the world of orthodontics. By developing products that correspond with this goal, LeoneAmerica hopes to meet the needs of the orthodontists and pediatric dentists working to optimize treatment for their patients. Wherever the OK Orthodontics for Kids logo is found, clinicians can rest assured that the particular product is appropriate and recommended especially for pediatric treatment. To get a copy of the new OK Orthodontics for Kids catalog, call us at (805) 487-9860.



**PEDODONTIC BANDS**

Designed to respond to the current needs of pediatric orthodontics and allow for early treatment of patient with mixed or deciduous teeth. Made of biomedical stainless steel in a softer medium temper, the pedodontic bands are designed for the particular anatomy and shape of the deciduous molars. These bands feature an easy fit to the pyramid shape of the primary molars, and are useful in the construction of palatal expanders or space maintainers. The laser etched identification number makes distinguishing between the ten sizes of upper or lower bands a cinch.



<b>E6100-00</b>	Upper - PU sizes 1-10	5
<b>E6500-00</b>	Lower - PL sizes 1-10	5

**E6165-96  
ASSORTMENT PEDODONTIC BANDS UNIVERSAL - KIT**

This kit is composed of 5 universal pedodontic bands per size for both maxillary and mandibular for a total of 100 pieces. The tray is not autoclavable.





**NEW**



The **MEMORIA®** Leaf Spring and Leaf Self Expanders are an evolution in the design of previous spring-loaded expanders. These innovative solutions rely on the flexible properties of nickel titanium springs to release calibrated and constant forces throughout treatment, providing expansion without the need for any patient compliance. The standard Leaf is activated periodically by the clinician in order to reload the springs, while the Leaf Self requires no re-activation at all. Both members of the series are available in either 6 or 9mm, with forces of either 450 or 900g.

**LEAF** *self* **EXPANDER®**



	11 mm	4 mm	∅	↔	↔	↔
			arms	body		
	6	450				
<b>A2705-06</b>	4 springs	450 g approx	1,5 mm	12 mm	6 mm	
	6	900				
<b>A2706-06</b>	4 springs	900 g approx				
	9	450				
<b>A2705-09</b>	6 springs	450 g approx	1,5 mm	12 mm	9 mm	
	9	900				
<b>A2706-09</b>	6 springs	900 g approx				

**LEAF** **EXPANDER®**



	11 mm	4 mm	∅	↔	↔	↻	activation turns
			arms	body			for maximum expansion limit
	6	450					
<b>A2703-06</b>	2 springs	450 g approx	1,5 mm	12 mm	6 mm	0,4 mm	30
	6	900					
<b>A2704-06</b>	2 springs	900 g approx					
	9	450					
<b>A2703-09</b>	3 springs	450 g approx	1,5 mm	12 mm	9 mm	0,4 mm	45
	9	900					
<b>A2704-09</b>	3 springs	900 g approx					

## CLASS II CORRECTOR (GOTO1®)

The non-compliance device for Class II correction. The small sizes of the device all for optimal patient comfort while the constant and light force delivered by the spring MEMORIA 200g, located inside the plunger, stimulates the mandibular advancement. In this new version, the fluidity of movement and the strength of the telescopic mechanism are improved. The packages include all the needed parts for the application of a bilateral correction device.

Available in three lengths with both mesial and distal fitting to the upper molar tube, for a total of six possible positions in the mouth. When used in the mesial position, it can be placed on a single direct-bonding tube, without the need for a band of round tube.




**mesial support**

to vestibular tube

**distal anchorage**


to vestibular tube



Two possible positions

1 with anterior eyelet

2 with anterior hook



Available in two versions

WITH ANTERIOR EYELET	WITH ANTERIOR HOOK	Type	Length at full expansion (passive configuration)	Length at full compression (including 2 mm of additional compensation)	Maximum safety stroke of spring	Maximum opening allowed to prevent disassembly
M2302-00	M2301-00	EXTRA-SHORT	31 mm	25 mm	7,5 mm	48 mm*
M2302-01	M2301-01	SHORT	35 mm	28 mm	10 mm	56 mm*
M2302-02	M2301-02	LONG	40 mm	33 mm	10 mm	62 mm*

\*Distance between the 2 eyelets or hooks of anchorage

# WHERE TO FIND US



## WWW.LEONEAMERICA.COM

LeoneAmerica Dental Products, Inc is delighted to introduce our upgraded website, [www.leoneamerica.com](http://www.leoneamerica.com). Products can be searched by orthodontics, pediatrics, or sleep apnea specialties, making the task of finding the perfect treatment solution even easier. We are so pleased to introduce this site as an interactive and useful tool for our clinical and laboratory customers.



## THROUGH TRUSTED SELLERS

Leone products are available through a number of trusted dealers throughout the United States. For more information on how to purchase Leone products, please call us at (805) 487-9860 or visit us online at [www.leoneamerica.com](http://www.leoneamerica.com)



## LeoLab USA

Certified LeoLab partners work closely with LeoneAmerica to bring you the best of research, innovation, and quality. The collaboration between specialized orthodontic experience and top-tier orthodontic manufacturing help to bring together some of the most precise orthodontic appliances on the market today.



## ASSOCIATION OF ORTHODONTIC LABORATORY PROFESSIONALS

LeoneAmerica Dental Products, Inc is proud to announce itself as a corporate sponsor of the Association of Orthodontic Laboratory Professionals. This collaborative group focuses on fostering relationships between orthodontic professionals, and moving the field of orthodontics together as a whole.

## ASK FOR LEONEAMERICA NEW CATALOG 2019

### Contact Info:

(805) 487-9860 | (800) 235-4639 | FAX (805) 834-0124  
[info@leoneamerica.com](mailto:info@leoneamerica.com) | [www.leoneamerica.com](http://www.leoneamerica.com)  
1250 Stellar Dr. Oxnard, CA 93033-2404

LeoneAmerica is the exclusive distributor of Leone Orthodontic products in the United States, Canada, and Puerto Rico.



LEONEAMERICA  
**ORTHO NEWS**

A PERIODICAL FOR ORTHODONTIC PROFESSIONALS

**SPRING/SUMMER 2019**  
VOL 1 #41

PRSR STD  
MAIL  
U.S. POSTAGE  
**PAID**  
OXNARD, CA  
PERMIT NO. 539

Booth #801

**May 3-7**



**American  
Association of  
Orthodontists®**

**2019**

**Los Angeles**



**MAY 3-7, 2019**

LOS ANGELES  
CONVENTION CENTER  
LOS ANGELES, CA

Booth #103

**AADSM  
ANNUAL  
MEETING**

**SAN ANTONIO: JUNE 7-9**

**19**



**JUNE 7-9, 2019**

AADSM ANNUAL SESSION  
MARRIOTT RIVERCENTER  
SAN ANTONIO TX

Booth #1717



**MAY 23-25, 2019**

AAPD ANNUAL SESSION  
HILTON CHICAGO  
CHICAGO, IL