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## PROSTHETIC TORQUE WRENCH 25 Ncm – INSTRUCTIONS FOR USE

**Warnings and restrictions.** The herewith furnished information is NOT intended for the patient, but exclusively for the medical professional, who will have to evaluate the contents in reference to every specific case to be treated, in order to furnish suitable information to the patient about the benefits and eventual side effects or complications. The prosthetic torque wrench must only be used by Dentists or Dental Surgeons experienced in implant techniques as well as in the related methodologies and manufacture of prosthesis on implants. The hereunder instructions for use are finalized to underline the risks, the contraindications and the collateral effects of general character (a reference to all the possible circumstances of use is not possible) and well known at the moment of the placing on the market of the device. Nevertheless the use of the prosthetic torque wrench requires a suitable knowledge of the prosthetic procedures illustrated on the Catalogue of the XCN® Implant System and on the web site [www.leone.it](http://www.leone.it). Considering that errors and deficiencies in patient evaluation, preoperative diagnosis, in the planning and the carrying out of the treatment may lead to implant loss, before use, it is strongly recommended to attend specific training courses with the purpose to reach a high degree of propaedeutical knowledge and practice in the use. It is not ascribed to Leone the duty of monitoring the procedures of use of the products for the implantprosthetic treatment; such procedures must integrally be reverted under the responsibility of the clinician. Leone takes no responsibility for the possible misuse of the device even neither for failure, adverse reactions or damages occurred to the patient or to the dentist following the use of the product. In any case the use of exclusively original components, accessories and instruments in conformity to the instructions is strongly recommended. **Description of the product.** The torque wrench is made of stainless steel. It is delivered together with an instrument lubricant. **Intended use.** The torque wrench is used to tighten connecting screws and the Conic adapter of the XCN® Implant System with a torque of 25 Ncm. **Utilization.** Detailed information regarding the use of the prosthetic torque wrench are available in the above mentioned Prosthetic Procedures. The torque wrench is supplied non-sterile, can be disassembled and must be sterilized by autoclave before use according to this instructions for use. The sterilization has to follow a validate procedure which is constantly monitored and follows ISO 17665-1 "Sterilization of health care products – moist heat – Part 1: requirements for the development, validation and routine control of a sterilization process for medical devices". It is recommended to sterilize the instruments and accessories with steam autoclave using the following parameters: temperature: 121°C (250°F), pressure: 1 atm, minimum exposure time: 20 minutes. The sterilization packaging must be appropriate for the device and conforming to EN 868-5 "Packaging for terminally sterilized medical devices - Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods".

We advise against the use of a hot air steriliser because it can lead to ageing of the spring and subsequently bring about a change of the torque value.

Before each use the device must be checked upon it's functioning.

The two black lines on the instrument have to be aligned to be in working position (Fig. A). No space must be visible between the two lines (Fig. B). If necessary, turn the handle clockwise to reach the right position.

Fig. A: correct position


Fig. B: wrong position

The word 'IN' on the cover of the torque wrench shows the position of the wrench that is used for tightening; by turning the device around, the word 'OUT' is used for loosening screws. **Recommendation.** The persons in charge for the use and maintenance of this dental instrument should monitor any deterioration of the tightening and torque mechanism of the device. **During assembly, it is essential not to mix the various components belonging to different torque wrenches because the components are not interchangeable, as each instrument is calibrated individually.** Components cannot be sold separately.

Fig. C: Individual parts of the torque wrench

- 1 ratchet wheel
- 2 ratchet
- 3 cover
- 4 ratchet head
- 5 cover screw (fixed within the ratchet head)
- 6 spring with guide pin
- 7 handle with hexagon head

### Cleaning the torque wrench

When used in situations that do lead to organic residues (blood, secretions, tissue remnants), the torque wrench must be disassembled completely (Fig. D to F) and placed in a suitable bath of disinfection in accordance with the recommendations of the manufacturer. This operation facilitates cleaning because dry residues cause corrosion. After cleaning, thoroughly rinse the parts with water and use a soft nylon brush to rub internal and external surfaces of the various parts of the torque wrench. During the cleaning process,

avoid all contact between each part of the torque wrench. Before sterilisation, the torque wrench must be completely assembled (Fig. G to K).

### Disassembling the prosthetic torque wrench

Fig. D: Unscrew completely the handle (7) and extract the spring with the guide pin (6). If necessary, use (gently) the hexagon head of the handle to extract the spring. Do not separate the guide pin from the spring (6).

Fig. D

Fig. E: Unscrew the cover screw (5) with the hexagon head of the handle. The screw is fixed within the ratchet head

Fig. E


Fig. F: Exert a light pressure in order to remove the cover (3) from the ratchet head (4). Take out the ratchet wheel (1) and the ratchet (2).

Fig. F

## Assembling the prosthetic torque wrench

Lubrication:

After cleaning, dry the single parts, assemble the instrument and moderately lubricate the functional contact areas as shown in Fig. G to J.

Remove any trace of the lubricant from the outer surface of the torque wrench. Use only the instrument lubricant delivered with the torque wrench.

Fig. G, H: Insert first the ratchet wheel (1) and then the ratcheted (2) into the ratchet head.


Fig. G

Fig. H

Fig. I: Put the cover (3) on the ratchet head and turn the torque wrench around while holding the cover. Fix the cover screw (5) tightly with the hexagon head of the handle.

Fig. I

Fig. J, K: Insert spring with guide pin into the ratchet head and screw the handle onto the ratchet head all the way to the stop. The two black lines on the instrument have to be aligned (Fig. K). Check the function by activating the ratchet wheel. If there is a regular ratchet noise, the instrument operates correctly.


Fig. J

Fig. K

**Contraindications. Side effects. Complications and interactions.** All the pathological conditions, both local and systemic, which contraindicate an intervention of oral surgery, have to be considered contraindications valid for instruments and accessories used for the application of dental implants and restorations as well. In the evaluation of the patient and the use of the product, Dentists and Dental Surgeons must contemplate contraindications, side effects, complications and interactions that can be generally applicable with oral surgery.

For further details regarding cleaning and sterilization please refer to the “Instruction for instruments sterilization” on the Leone website [www.leone.it](http://www.leone.it) section Service/Quality.

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Below are the symbols used, where applicable, on the label of surgical instruments and accessories.