

D.B. BRACKETS F1000PLUS



**F1000PLUS
ROTH SYSTEM**

						
			torque	ang.	.022"x.030"	
			+12°	+5°		F1000-11
						F1000-21
			+8°	+9°		F1000-12
						F1000-22
			0°	+11°		F1002-13
						F1002-23
			-7°	0°		F1002-14
						F1002-24
			-7°	0°		F1002-15
						F1002-25
			-2°	0°		F1002-41
						F1002-31
			-2°	0°		F1002-42
						F1002-32
			-11°	+5°		F1002-43
						F1002-33
			-17°	0°		F1002-44
						F1002-34
			-22°	0°		F1002-45
						F1002-35



**F1000PLUS
MBT* SYSTEM**

						
			torque	ang.	.022"x.030"	
			+17°	+4°		F1001-11
						F1001-21
			+8°	+9°		F1001-12
						F1001-22
			-8°	+8°		F1001-13
						F1001-23
			-7°	0°		F1000-14
						F1000-24
			-7°	0°		F1000-15
						F1000-25
			-6°	0°		F1001-41
						F1001-31
			-6°	0°		F1001-42
						F1001-32
			-6°	+3°		F1001-43
						F1001-33
			-12°	0°		F1000-44
						F1000-34
			-17°	0°		F1000-45
						F1000-35



**F1000PLUS
ROTH SYSTEM**

	20 brackets 1 case	200 brackets 10 cases
.022"	F1000-93	F1001-93



**F1000PLUS
MBT* SYSTEM**

	20 brackets 1 case	200 brackets 10 cases
.022"	F1000-92	F1001-92

*MBT is a Trademark of 3M Unitek.
The orthodontic brackets illustrated in this catalogue are not intended to be a duplication of any other existing system nor does Leone SpA imply that they are endorsed by the above mentioned doctors or Schools in any form.

D.B. BRACKETS **F1000PLUS**



F1000PLUS
DAMON SYSTEM

						
torque		ang.			.022"x.030"	
			+12°	+5°	1 1	F1000-11 F1000-21
			+8°	+9°	2 2	F1000-12 F1000-22
			+7°	+8°	3 3	F1000-13 F1000-23
			-7°	0°	4 4	F1000-14 F1000-24
			-7°	0°	5 5	F1000-15 F1000-25
			-2°	0°	1 1	F1000-41 F1000-31
			-2°	0°	2 2	F1000-42 F1000-32
			+7°	+3°	3 3	F1000-43 F1000-33
			-12°	0°	4 4	F1000-44 F1000-34
			-17°	0°	5 5	F1000-45 F1000-35



F1000PLUS
DAMON SYSTEM

	20 brackets 1 case	200 brackets 10 cases
.022"	F1000-91	F1001-91



**INTERACTIVE SL
D.B. BRACKETS**

Pack of 5



**INTERACTIVE SL
ROTH SYSTEM**

						.022"x.030"
			+12°	+5°	$\frac{1}{1}$	F1100-11 F1100-21
			+8°	+9°	$\frac{2}{2}$	F1100-12 F1100-22
			-2°	+13°	$\frac{3}{3}$	F1100-13 F1100-23
			-7°	0°	$\frac{4}{4}$	F1100-14 F1100-24
			-7°	0°	$\frac{5}{5}$	F1100-15 F1100-25
			0°		$\frac{1}{1}$	F1100-41 F1100-31
			0°		$\frac{2}{2}$	F1100-42 F1100-32
			-11°	+7°	$\frac{3}{3}$	F1100-43 F1100-33
			-17°	0°	$\frac{4}{4}$	F1100-44 F1100-34
			-22°	0°	$\frac{5}{5}$	F1100-45 F1100-35



**INTERACTIVE SL
MBT* SYSTEM**

						.022"x.030"
			+17°	+4°	$\frac{1}{1}$	F1102-11 F1102-21
			+10°	+8°	$\frac{2}{2}$	F1102-12 F1102-22
			-7°	+8°	$\frac{3}{3}$	F1102-13 F1102-23
			-7°	0°	$\frac{4}{4}$	F1100-14 F1100-24
			-7°	0°	$\frac{5}{5}$	F1100-15 F1100-25
			-6°	0°	$\frac{1}{1}$	F1102-41 F1102-31
			-6°	0°	$\frac{2}{2}$	F1102-42 F1102-32
			-6°	+3°	$\frac{3}{3}$	F1102-43 F1102-33
			-12°	+2°	$\frac{4}{4}$	F1102-44 F1102-34
			-17°	+2°	$\frac{5}{5}$	F1102-45 F1102-35



**INTERACTIVE SL
ROTH SYSTEM**

	20 brackets 1 case	200 brackets 10 cases
.022"	F1100-91	F1101-91



**INTERACTIVE SL
MBT* SYSTEM**

	20 brackets 1 case	200 brackets 10 cases
.022"	F1100-92	F1101-92

*MBT is a Trademark of 3M Unitek.
The orthodontic brackets illustrated in this catalogue are not intended to be a duplication of any other existing system nor does Leone SpA imply that they are endorsed by the above mentioned doctors or Schools in any form.